

CONCUSSION CARE - NEXT STEPS

Using the ImPACT Passport® App

You have recently completed **ImPACT, ImPACT Pediatric, or ImPACT Quick Test** as part of a comprehensive concussion care evaluation.

Now what?

Your test scores are saved in a clinical report intended for a trained healthcare provider. It's important that **concussion care providers can access your clinical report when you need help** as you won't have a direct access to it due to its clinical nature. Your ImPACT Passport ID is a key that helps unlock your clinical data by a trained healthcare provider.

For that reason, ImPACT Applications developed **ImPACT Passport - your personal concussion care app.**



ImPACT Passport allows you to:

- Store your unique testing ID
- Learn more about signs and symptoms of a concussion
- Record symptoms
- Grant healthcare providers access to your data




CONCUSSION CARE - NEXT STEPS

Using the ImPACT Passport® App

Where to get your ImPACT Passport ID?

Upon ImPACT completion, you'll get an ImPACT Passport ID, which you can choose to send to your preferred email. Make sure to enter your ImPACT Passport ID into your ImPACT Passport app.



ImPACT Passport ID: W3KU-2IV3-MLTX

Upon ImPACT Pediatric or ImPACT Quick Test completion, ask the healthcare provider who administered the test to give you your ImPACT Passport ID. Make sure to enter your ID into your ImPACT Passport app.

How to get ImPACT Passport?




ImPACT Passport is available for download at Apple App Store and Google Play


What to do with ImPACT Passport?




Save your ImPACT Passport ID to grant access to your clinical report



Learn about concussion signs and symptoms



Report symptoms so your healthcare provider can review them real-time



Find a concussion care provider near you

