



The most well-researched, computerized concussion assessment tool just got even better with the release of new product developments.

ImPACT Version 4

What is ImPACT Version 4?

ImPACT Version 4 continues to use the same clinically proven tests as in the previous versions and includes new features, updated normative data, and expanded use cases.

What's new in ImPACT Version 4?

Broader application to concussion care:

- > Applicable age range for ImPACT extended to 80 (the new age range is 12-80)
- > Test takers given the option to use a trackpad or mouse

Updated normative values:

- > Updated normative values for ages 12-80
- > Separate normative values for tests taken on trackpad or mouse
- > Revised invalidity indicators to better reflect the new normative data

Streamlined clinical report:

- > Addition of the **Two-Factor score** (memory and speed)
- > Removal of Cognitive Efficiency Index (CEI)

How do I get started with ImPACT Version 4?

ImPACT Applications customers can immediately access the new features available in ImPACT Version 4 at no additional cost.

Have questions or need help getting started? Check out the [ImPACT Version 4 FAQs](#) or email us at info@impacttest.com.